

**Digital dry bath LLG-uniBLOCKTHERM****LLG Labware****Digital dry bath LLG-uniBLOCKTHERM**

Suitable for use in microbiological, histological, clinical, industrial and environmental laboratories. Various types of optional heating blocks are available as accessories for the LLG-uniBLOCKTHERM.

- Digital control with LED display
- Can be used up to 120 °C
- For timer and continuous operation
- Can be used with different heating blocks
- 1 x block unit
- With overheating protection

Scope of supply: LLG-uniBLOCKTHERM, removal tool, lid

Please order heating blocks separately.

Specifications

Temperature range: RT ... 120 °C

Temperature uniformity: ± 0.5 °C

Temperature accuracy: ± 0.5 °C

Timer: 0 ... 99:59:00 min

Block material: Aluminium

Dimensions (W x D x H): 175 x 290 x 85 mm

Weight: 3 kg

Power supply: 220 ... 240 V, 50/60 Hz

Warranty: 3 years



Code	Description	Packaging
LLG06263460	LLG-uniBLOCKTHERM, Block thermostat	1 pz.
LLG06263461	LLG-uniBLOCKTHERM, UK-Plug	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

