



Bottletop dispenser LLG-uniTOPDISPENS 2

LLG Labware

Bottletop dispenser LLG-uniTOPDISPENS 2

The combination of a springless PTFE valve and **PTFE piston** with chemically resistant ETP O-ring enables dispensing of both organic reagents as well as bases, thus, eliminating the need to buy different dispensers for different applications. The improved design of the uniTOPDISPENS 2 bottle top dispenser including a **re-circulation valve** allows a smooth, precise and reliable bubble free dispensing without any loss of reagent. The specially designed bottle adaptor which permits a 360° rotation of the unit, the easy volume setting as well as the flexible delivery nozzle underline the user friendliness of this bottletop dispenser. Available in six volume ranges.

- Springless PTFE valve, PTFE piston with chemically resistant ETP O-ring enables **high chemical compatibility**
- Glass cylinder from borosilicate glass 3.3
- Re-circulation valve enables easy priming and bubble free dispensing
- Smooth, reproducible and effortless dispensing
- 360° rotation possible due to specially designed bottle adapter
- Flexible and adjustable delivery nozzle provides **dispensing angle flexibility** to the user
- Specially designed volume adjustment knob for quick and easy volume setting
- The telescopic extractable intake tube can be easily adjusted to varying bottle heights**
- Nozzle cap reduces contact with air and prevents drying and dripping if the tip cover is attached
- Five bottle adaptors included (28 mm, 32 mm, 37 mm, 40 mm and 45 mm). Standard thread: 30 mm.
- Fully autoclavable at 121 °C, 15 psi for 10 to 15 min
- Individually calibrated in accordance with ISO 8655 standards and each unit is supplied with an individual calibration certificate
- You can find the appropriate compatibility list on www.llg-labware.com in the Downloads section.



-3 years warranty

Code	Description	Packaging
LLG06282275	Bottletop dispenser uniTOPDISPENS 2, uniTOPDISPENS 2, 0.25 - 2.5 ml, Graduations 0.05 ml	1 pz.
LLG06282276	Bottletop dispenser uniTOPDISPENS 2, uniTOPDISPENS 2, 0.5 - 5 ml, Graduations 0.1 ml	1 pz.
LLG06282278	Bottletop dispenser uniTOPDISPENS 2, uniTOPDISPENS 2, 2.5 - 30 ml, Graduations 0.5 ml	1 pz.
LLG06282279	Bottletop dispenser uniTOPDISPENS 2, uniTOPDISPENS 2, 5 - 60 ml, Graduations 1.0 ml	1 pz.



Code	Description	Packaging
LLG06282277	Bottletop dispenser uniTOPDISPENS 2, uniTOPDISPENS 2, 1 - 10 ml, Graduations 0.2 ml	1 pz.
LLG06282280	Bottletop dispenser uniTOPDISPENS 2, uniTOPDISPENS 2, 10 - 100 ml, Graduations 2.0 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

