

**LLG Syringe Filter CA, Cellulose acetate****LLG Labware****LLG Syringe Filter CA, Cellulose acetate**

We offer a line of syringe filters especially designed to provide efficient filtration of a broad variety of fluids, solvents, aqueous or inorganic solutions. Our syringe filters cover most applications in HPLC, pharmaceutical, environmental, biotechnology and food and beverage testing laboratories.

- Housing injected in Polypropylene or MABS
- Multifunctional Syringe filters equipped with Luer-Lock connections on both sides or female Luer-Lock input and male Luer-Slip output
- Suitable for all syringes with Luer connection
- Pre-sterilised or non-sterile
- Sterile products in individual hard blister
- Superior pressure stability

Hydrophilic membrane for filtration of aqueous solutions.

This membrane features excellent shape stability in aqueous solutions and a very low binding capacity for proteins (21µg per 25 mm Filter). Ideal for use with biological macromolecules, water soluble oligomers and polymers. With Luer-Lock input and Luer-Lock or Luer-Slip output.



Code	Description	Packaging
LLG06285694	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 13 mm, Pore size 0.8 µm, Housing Acrylic, green	50 pz.
LLG09055504	LLG Syringe Filter CA, Cellulose acetate, Type non sterile, Filter diam. 13 mm, Pore size 0.8 µm, Housing Acrylic, green	500 pz
LLG06285704	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 25 mm, Pore size 0.45 µm, Housing Acrylic, transparent*	50 pz.
LLG06285705	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 25 mm, Pore size 0.8 µm, Housing Acrylic, transparent*	50 pz.
LLG09055500	LLG Syringe Filter CA, Cellulose acetate, Type non sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing Acrylic, blue	500 pz
LLG09055502	LLG Syringe Filter CA, Cellulose acetate, Type non sterile, Filter diam. 13 mm, Pore size 0.45 µm, Housing Acrylic, yellow	500 pz
LLG09055510	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 13 mm, Pore size 0.2 µm, Housing Acrylic, blue	50 pz
LLG09055512	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 13 mm, Pore size 0.45 µm, Housing Acrylic, yellow	50 pz
LLG06285699	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 25 mm, Pore size 0.8 µm, Housing Acrylic, green	50 pz.
LLG06285703	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 25 mm, Pore size 0.2 µm, Housing Acrylic, transparent*	50 pz.
LLG07970389	LLG Syringe Filter CA, Cellulose acetate, Type non sterile, Filter diam. 25 mm, Pore size 0.8 µm, Housing Acrylic, green	500 pz
LLG09055501	LLG Syringe Filter CA, Cellulose acetate, Type non sterile, Filter diam. 25 mm, Pore size 0.2 µm, Housing Acrylic, blue	500 pz
LLG09055503	LLG Syringe Filter CA, Cellulose acetate, Type non sterile, Filter diam. 25 mm, Pore size 0.45 µm, Housing Acrylic, yellow	500 pz
LLG09055511	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 25 mm, Pore size 0.2 µm, Housing Acrylic, blue	50 pz



Code	Description	Packaging
LLG09055513	LLG Syringe Filter CA, Cellulose acetate, Type sterile, Filter diam. 25 mm, Pore size 0.45 µm, Housing Acrylic, yellow	50 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

