



Isolating box with lid, Neopor

Isolating box with lid

The insulation boxes are suitable for the transport and storage of temperature-sensitive products.

- To maintain cold chains
- Stackable

Neopor® is characterised by its high hot and cold thermal insulation properties, compressive strength, shock absorbing properties, low weight and insensitivity to moisture. In addition, it has infrared absorbers and reflectors that reduce thermal conductivity and give the box its silver-grey colour.



Code	Description	Packaging
LLG07656681	Isolating box with lid, Isolating baffle ring, 6.2 L	1 pz
LLG07656680	Isolating box with lid, Isolating box with lid, 12.5 L	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

