

**Tanks Optima series, stainless steel****Grant****Tanks Optimaä series, stainless steel**

High efficiency stainless steel (ST) tanks (range -15 °C to 200 °C). Drain tap included on ST12, ST18, ST26, ST38. Raised feet for carrying repositioning and retort access. Accessory polypropylene spheres or a lid must be used at temperatures between 60 °C and 100 °C. Above 100 °C only a lid must be used. Accessory cooling is required for temperatures below ambient.

Please order lid and PP spheres separately.



Code	Description	Packaging
LLG09905722	Tanks Optima series, stainless steel, Type ST38, 38 l, External dim. (W x D x H) 330 x 730 x 260 mm, Internal dim. (W x D x H) 300 x 690 x 200 mm	1 pz.
LLG09905720	Tanks Optima series, stainless steel, Type ST18, 18 l, External dim. (W x D x H) 330 x 540 x 215 mm, Internal dim. (W x D x H) 300 x 505 x 150 mm	1 pz.
LLG09905721	Tanks Optima series, stainless steel, Type ST26, 26 l, External dim. (W x D x H) 330 x 540 x 265 mm, Internal dim. (W x D x H) 300 x 505 x 200 mm	1 pz.
LLG09905718	Tanks Optima series, stainless steel, Type ST5, 5 l, External dim. (W x D x H) 180 x 330 x 200 mm, Internal dim. (W x D x H) 150 x 300 x 150 mm	1 pz.
LLG09905719	Tanks Optima series, stainless steel, Type ST12, 12 l, External dim. (W x D x H) 330 x 360 x 215 mm, Internal dim. (W x D x H) 300 x 325 x 150 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

