



Test tubes, Soda-lime glass, heavy wall

Scherf

Test tubes, Soda-lime glass, heavy wall

Rimless, with round bottom (0.7 to 1.0 mm wall thickness).



Code	Description	Packaging
LLG06258821	Test tubes, Soda-lime glass, heavy wall, Diam. 12 mm	250 pz.
LLG06258822	Test tubes, Soda-lime glass, heavy wall, Diam. 12 mm	250 pz.
LLG06258823	Test tubes, Soda-lime glass, heavy wall, Diam. 14 mm	250 pz.
LLG06258824	Test tubes, Soda-lime glass, heavy wall, Diam. 14 mm	250 pz.
LLG06258825	Test tubes, Soda-lime glass, heavy wall, Diam. 16 mm	250 pz.
LLG06258826	Test tubes, Soda-lime glass, heavy wall, Diam. 16 mm	250 pz.
LLG06258827	Test tubes, Soda-lime glass, heavy wall, Diam. 18 mm	250 pz.
LLG09190200	Test tubes, Soda-lime glass, heavy wall, Diam. 8 mm	100 pz.
LLG09190202	Test tubes, Soda-lime glass, heavy wall, Diam. 10 mm	100 pz.
LLG09190205	Test tubes, Soda-lime glass, heavy wall, Diam. 10 mm	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

