



Rotator Genie, with platform

Rotator Genieä, with platform

Holds a wide variety of tubes, flasks, microplates, and odd-shaped vessels. Available in two interchangeable platform configurations. Magnetic platform features unique spring clip plate system. The standard platform is designed to turn 90 degrees allowing different mixing motions. Designed for cold room and incubator use. Will easily fit in most incubators. Uniquely designed for a wide selection of motions, which can be adjusted from gentle rolling to forceful end over end mixing. Continuous-duty quality motor provides smooth, quiet operation.

Specifications

Speed: 5rpm to 35rpm

Environmental: ATEX-Installation Category 2, 0°C to 38°C, 95% max. relative humidity

Overall dimensions (W x D x H): 400 x 175 (290 platform horizontal) x 400mm

Dimensions magnetic platform: 290 x 315mm

Dimensions standard platform: 260 x 260mm

Weight: 5.2kg



Code	Description	Packaging
LLG09730102	Rotator Genie, magn.platform 230 V, UK-plug	1 pz.
LLG09730103	Rotator Genie, magn.platform 230 V, swiss plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

