



## Rotator Genie, with platform

### Rotator Genieä, with platform

Holds a wide variety of tubes, flasks, microplates, and odd-shaped vessels. Available in two interchangeable platform configurations. Magnetic platform features unique spring clip plate system. The standard platform is designed to turn 90 degrees allowing different mixing motions. Designed for cold room and incubator use. Will easily fit in most incubators. Uniquely designed for a wide selection of motions, which can be adjusted from gentle rolling to forceful end over end mixing. Continuous-duty quality motor provides smooth, quiet operation.

#### Specifications

**Speed: 5rpm to 35rpm**

**Environmental: ATEX-Installation Category 2, 0°C to 38°C, 95% max. relative humidity**

**Overall dimensions (W x D x H): 400 x 175 (290 platform horizontal) x 400mm**

**Dimensions magnetic platform: 290 x 315mm**

**Dimensions standard platform: 260 x 260mm**

**Weight: 5.2kg**



| Code        | Description                                 | Packaging |
|-------------|---|-----------|
| LLG09730102 | Rotator Genie, magn.platform 230 V, UK-plug | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

