

**Jerrycan with integrated spout, HDPE****Hünersdorff****Jerrycan with integrated spout, HDPE**

Jerrycan in stable design made of food safe polyethylene (HDPE)

- UV protection for long life
- Especially suitable for drinking water and food
- Space-saving and stackable
- Integrated spout with wide thread
- Fill opening 65mm diameter
- With threaded plug and wide thread 22mm



Code	Description	Packaging
LLG09002611	Jerrycan with integrated spout, HDPE, 10 Litres, Width 350 mm, Depth 165 mm, Height 310 mm	1 pz.
LLG09002612	Jerrycan with integrated spout, HDPE, 20 Litres, Width 350 mm, Depth 165 mm, Height 495 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

