

**Universal Ovens UN and UF****Memmert****Universal Ovens UN and UF**

The all-round genius among heating ovens covers a multitude of applications, ideally at temperatures above +50 °C. The fresh air is preheated in the pre-heating chamber to avoid temperature fluctuations.

SingleDISPLAY: ControlCOCKPIT with one colour display

You can find a detailed description of the "SingleDISPLAY" model at www.memmert.com



Code	Description	Packaging
LLG06266314	Universal cabinet UF1060 +20...+300°C, 1060 ltr. laboured air circulation	1 pz.
LLG09869681	Universal cabinet UN30 +20...+300°C, 32 ltr. natural air circulation	1 pz.
LLG09869684	Universal cabinet UN110 +20...+300°C, 108 ltr. natural air circulation	1 pz.
LLG09869702	Universal cabinet UF260 +20...+300°C, 256 ltr. laboured air circulation	1 pz.
LLG09869703	Universal cabinet UF450 +20...+300°C, 449 ltr. laboured air circulation	1 pz.
LLG09869704	Universal cabinet UF750 +20...+300°C, 749 ltr. laboured air circulation	1 pz.
LLG09869701	Universal cabinet UF160 +20...+300°C, 161 ltr. laboured air circulation	1 pz.
LLG09869682	Universal cabinet UN55 +20...+300°C, 53 ltr. natural air circulation	1 pz.
LLG09869683	Universal cabinet UN75 +20...+300°C, 74 ltr. natural air circulation	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

