



## Hand Disinfectant Softa-Man

### Hand Disinfectant Softa-Man®

Ready to use alcohol-based solution for gentle hand disinfection.

- Effective against bacteria including fungi, lipophile viruses including HBV and HIV
- Gentle cleansing agents and non-greasy emollients make it non-allergenic
- Dermatologically approved
- DGHM- and BGA-listed
- For normal to dry skin

This product is a disinfectant containing biocides. Please inquire its availability and legal information.



Code	Description	Packaging
LLG09191213	Hand Disinfectant Softa-Man®, Type Softa-Man® Dispenser, 1000 ml	1 pz.
LLG09191216	Hand Disinfectant Softa-Man®, Type Softa-Man® ViscoRub Dispenser, 1000 ml	1 pz.
LLG09191215	Hand Disinfectant Softa-Man®, Type Softa-Man® Canister, 5000 ml	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

