

**LLG-Laboratory jack, stainless steel****LLG Labware****LLG-Laboratory jack, stainless steel**

- Plates made of stainless steel.
- Top plate with Ø 10 mm hole.
- With non-skid rubber feet.
- Acc. to DIN 12897.



Code	Description	Packaging
LLG09117220	LLG-Laboratory jack, stainless steel	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

