



Cells for Fluorescence Measurement, UV-range

Cells for Fluorescence Measurement, UV-range

Wavelength Range: 200nm - 2500nm

Optical material: Quartz glass High Performance



Code	Description	Packaging
LLG06051579	Cells for Fluorescence Measurement, UV-range, Type ultramicro, PE stopper, 3 window	1 pz.
LLG06081483	Cells for Fluorescence Measurement, UV-range, Type micro, PTFE lid	1 pz.
LLG06093874	Cells for Fluorescence Measurement, UV-range, Type semi-micro, PTFE lid	1 pz.
LLG06204497	Cells for Fluorescence Measurement, UV-range, Type macro, PTFE stopper, 4 window	1 pz.
LLG06302645	Cells for Fluorescence Measurement, UV-range, Type macro, clip ***	1 pz.
LLG06801662	Cells for Fluorescence Measurement, UV-range, Type macro, clip ***	1 pz.
LLG09144759	Micro cuvette 115F-QS, 10x2mm coat thic	1 pz.
LLG09144761	Cells for Fluorescence Measurement, UV-range, Type macro, PTFE lid, 4 window	1 pz.
LLG09144771	Cells for Fluorescence Measurement, UV-range, Type macro, PTFE stopper, 4 window	1 pz.
LLG09144781	Cells for Fluorescence Measurement, UV-range, Type semi-micro, PTFE lid	1 pz.
LLG09144791	Fluorescence cells, quark glasspatz len	1 pz.
LLG09144801	Cells for Fluorescence Measurement, UV-range, Type micro, PTFE lid	1 pz.
LLG07638526	Cells for Fluorescence Measurement, UV-range, Type macro, PTFE lid, 4 window	1 pz.
LLG09144390	Cells for Fluorescence Measurement, UV-range, Type ultramicro, PE stopper, 3 window	1 pz.
LLG09144391	Cells for Fluorescence Measurement, UV-range, Type ultramicro, PE stopper, 3 window	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

