

**Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps****Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps**

Fittings with integrated ferrule for SafetyCaps and SafetyWasteCaps.

Use the blind plugs to tighten connectors which are temporarily inactive or not in use.



Code	Description	Packaging
LLG04005445	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, Material PTFE, Colour white, Outer diam. 2.3 mm	10 pz.
LLG09042892	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, Material PFA, Colour blue, Outer diam. 3.2 mm	5 pz.
LLG04005447	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, PTFE, Colour white, Outer diam. 6.35 (1/4") mm	1 pz.
LLG06241792	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, Material PFA, Colour green, Outer diam. 1.6 mm	5 pz.
LLG09139890	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Blind plugs for capillary connections, PFA, Colour neutral, Outer diam. 1.6/2.3/3.2 mm	10 pz.
LLG06255961	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, Material PFA, Colour violet, Outer diam. 2.3 mm	5 pz.
LLG04005444	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, Material PTFE, Colour white, Outer diam. 1.6 mm	10 pz.
LLG04005446	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, Material PTFE, Colour white, Outer diam. 3.2 mm	10 pz.
LLG04005448	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Fittings for capillaries, PTFE, Colour white, Outer diam. 4.76 (3/16") mm	1 pz.
LLG04005880	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Blind plugs for capillary connections, PTFE, Colour white, Outer diam. 4.76 mm	10 pz.
LLG04005945	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Blind plugs for capillary connections, PTFE, Colour white, Outer diam. 4.76 mm	5 pz.



Code	Description	Packaging
LLG04005861	T PTFE pour tube diam ext:1,6/2,3/3,2 mm	
LLG04005883	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Blind plug for tubing connection, Material PTFE, Colour white, Outer diam. - mm	1 pz.
LLG04005504	Fittings and Ferrules for capillary connector for SafetyCaps / SafetyWasteCaps, Blind plugs for charcoal filter connections, PTFE/VITON, Colour white/black, Outer diam. - mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

