

**Collectors for tube connector for SafetyWasteCaps****S.C.A.T.****Collectors for tube connector for SafetyWasteCaps**

Collectors for additional parts to add more tubings and capillaries. Delivered with fittings.

Further configurations available on request.



Code	Description	Packaging
LLG04005865	Collectors for tube connector for SafetyWasteCaps, 2 x connectors (90°) for 3 tubings diam. 6.4-9.0mm i.d., Material PTFE/PP	1 pz.
LLG04005866	Collectors for tube connector for SafetyWasteCaps, 2 x connectors (straight) for 2 capillaries diam. 2.3/3.2mm o.d., Material PTFE/PFA	1 pz.
LLG04005859	Collectors for tube connector for SafetyWasteCaps, 3 x connectors (90°) for 2 capillaries diam. 2.3/3.2mm o.d. and 1 tubing diam. 6.4-9.0mm i.d., PTFE/PFA/PP	1 pz.
LLG04005860	Collectors for tube connector for SafetyWasteCaps, 3 x connectors (straight) for 3 capillaries diam. 2.3/3.2mm o.d., Material PTFE/PFA	1 pz.
LLG06262289	Collectors for tube connector for SafetyWasteCaps, 3 x connectors (straight) for 3 capillaries diam. 6.4mm i.d., Material PTFE/PP	1 pz.
LLG04005976	Collectors for tube connector for SafetyWasteCaps, 8 x connectors (straight) for 7 capillaries diam. 2.3/3.2mm o.d. and 1 tubing diam. 6.4mm i.d., PTFE/PFA/PP	1 pz.
LLG09139889	Collectors for tube connector for SafetyWasteCaps, 8 x connectors (straight) for 8 capillaries diam. 2.3/3.2mm o.d., Material PTFE/PFA	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

