

**InKjel P infrared digestion system with 10 freely configurable programmes****behr****InKjel P infrared digestion system with 10 freely configurable programmes**

Rapid digestion system with 25 freely configurable programmes for power and digestion time and direct heating of the samples through high-quality quartz radiators (1500W).

- Quartz infra-red radiators peak at 830°C within 1 min.
- Particularly uniform heating at all sample positions

All InKjel systems are fitted completely with multi-level consoles, fume extraction unit, sample rack and glass digestion vessels.

Configurations available:

- 6 sample vessels each 250ml
- 12 sample vessels each 100ml
- 12 sample vessels each 250ml
- 4 sample vessels each 500ml
- 4 sample vessels each 750ml



Code	Description	Packaging
LLG09843756	InKjel P infrared digestion system with 10 freely configurable programmes, Type InKjel 1210 P, For 12 sample vessels each 100ml	1 pz.
LLG09843758	InKjel P infrared digestion system with 10 freely configurable programmes, Type InKjel 450 P, For 4 sample vessels each 500ml	1 pz.
LLG09843759	InKjel P infrared digestion system with 10 freely configurable programmes, Type InKjel 475 P, For 4 sample vessels each 750ml	1 pz.
LLG09843760	InKjel P infrared digestion system with 10 freely configurable programmes, Type InKjel 625 P, For 6 sample vessels each 250ml	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

