

**Homogeniser T 10 basic ULTRA-TURRAX****IKA****Homogeniser T 10 basic ULTRA-TURRAX®**

Small dispersing instrument for volumes from 0.5ml to 100ml (H₂O). It offers a wide speed range from 8000rpm to 30000rpm that enables users to work at high circumferential speeds even with small rotor diameters. A broad choice of dispersing elements guarantees a wide range of applications.

- Quick-connect coupling for changing dispersing elements easily
- Ideal for manual operation due to its light weight and ergonomic design
- Extremely mobile due to direct line power (no transformer required)
- Stainless steel dispersing elements can be cleaned quickly and easily as they can be dismantled without tools
- Plastic disposable dispersing elements are available in two sizes, particularly suitable for PCR analysis
- High performance drive ensures immense speed stability
- Quiet operation

Dispersing heads must be ordered separately.

Specifications

Motor rating input / output: 125W / 75W

Working volume range (H₂O): 0.5ml to 100ml

Viscosity max.: 5000mPas

Speed range: 8000rpm to 30000rpm

Speed display: scale

Speed control: stepless

Noise without element: 65dB(A)

Extension arm diam. / length: 8mm / 130mm

Dimensions (W x D x H): 56 x 66 x 178mm

Weight: 0.5kg

Power supply: 230V, 50/60Hz

Protection class DIN EN 60529: IP 30



Code	Description	Packaging
LLG09570940	ULTRA-TURRAX® T 10 basic Disperser, 220...240 V, 50/60 Hz	1 pz.
LLG09570941	ULTRA-TURRAX® T 10 basic S 2 Disperser, with UK plug, 220...240 V, 50/60 Hz	1 pz.
LLG09570942	ULTRA-TURRAX® T 10 basic S 3 Disperser, with Swiss plug 220...240 V, 50/60 Hz	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

