

**Homogeniser T 18 digital ULTRA-TURRAX****IKA****Homogeniser T 18 digital ULTRA-TURRAX®**

Competitively priced dispersing instrument for volumes from 1ml to 1500ml (H₂O) with digital speed display. It offers a wide speed range from 500rpm to 25000rpm that enables users to work at high circumferential speeds even with small rotor diameters. A broad choice of dispersing elements guarantees a wide range of applications.

- Digital speed display
- Electronic speed control
- Electronic overload protection
- Stainless steel dispersing elements can be cleaned quickly and easily
- Plastic disposable dispersing elements are available in two sizes
- Error code display
- Quiet operation

Dispersing heads must be ordered separately.

Specifications

Motor rating input / output: 500W / 300W

Working volume range (H₂O): 1ml to 1500ml

Viscosity max.: 5000mPas

Speed range: 500rpm to 25000rpm

Speed display: LED

Speed control: stepless

Noise without element: 75dB(A)

Extension arm diam. / length: 13mm / 160mm

Dimensions (W x D x H): 87 x 106 x 271mm

Weight: 2.5kg

Power supply: 200-240 V, 50/60Hz

Protection class DIN EN 60529: IP 20



Code	Description	Packaging
LLG09570943	ULTRA-TURRAX® T 18 digital Disperser, 230 V, 50/60 Hz	1 pz.
LLG09570944	ULTRA-TURRAX® T 18 digital S 2 Disperser, with UK plug 230 V, 50/60 Hz	1 pz.
LLG09570945	ULTRA-TURRAX® T 18 digital S 3 Disperser, with Swiss plug 230 V, 50/60 Hz	1 pz.
LLG09570946	ULTRA-TURRAX® T 18 digital Package Disperser, incl. 9.570 842 230 V, 50/60 Hz	1 pz.
LLG09570947	ULTRA-TURRAX® T 18 digital Package S 2 Disperser, incl. 9.570 842 230 V, 50/60 Hz, with UK plug	1 pz.
LLG09570948	ULTRA-TURRAX® T 18 digital package S 3 Disperser, incl. 9.570 842 230 V, 50/60 Hz, with Swiss plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

