

**Multi-position magnetic stirrers RO 5/10/15 series****IKA****Multi-position magnetic stirrers RO 5/10/15 series**

Without heating. 5, 10 or 15-position digital magnetic stirrers, designed for synchronous stirring. The magnetic coil technology provides noiseless and consistent stirring on all positions. Speed remains constant, even when load changes. The closed and compact design allows easy cleaning and protects the equipment against the penetration of liquids. Additionally, a removable transparent anti-slip mat is included.

- Adjustable reverse operation
- Eco-mode for a low self-warming of the surface
- Wear free magnetic coils
- Error Code Display
- Easy operation with touch keypad
- Suitable for continuous operation
- Adjustable speed in increments of 10rpm

Specifications

Max. stirring quantity per stirring position (H₂O): 0.4l

Distance between stirring places: 90mm

Deviation for individual stirring positions: 0%

Speed range: 0 to 1200rpm

Speed display: LED line

Set-up plate material: stainless steel 1.4301

Protection class acc. to DIN EN 60529: IP 40

Power supply: 100-240V, 50/60Hz



Code	Description	Packaging
LLG09720222	Multi-position magnetic stirrer RO 10 digital, with 10 stirring places, w/o heating	1 pz.
LLG09720208	Multi position magnetic stirrer RO 5 digital, with 5 stirring places, w/o heating	1 pz.
LLG09720226	Multi-position magnetic stirrer RO 15 digital, with 15 stirring places, w/o heating	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

