

**Magnetic stirrers/hotplates, multi-position, RT 5/10/15****IKA****Magnetic stirrers/hotplates, multi-position, RT 5/10/15**

Digital magnetic hotplate stirrers with 5, 10 or 15 positions, designed for synchronous heating and stirring. The magnetic coil technology provides noiseless and consistent stirring on all positions. The surface temperature is infinitely adjustable up to 120°C, producing a maximum medium temperature of 70°C (depending on type of vessel). Speed remains constant, even when load changes.

- Homogeneous temperature distribution across heating plate
- Adjustable reverse operation
- Eco-mode for a low self-warming of the surface
- Wear free magnetic coils
- Error code display
- Easy operation with touch keypad
- Suitable for continuous operation
- Adjustable speed in increments of 10rpm

Specifications

Max. stirring quantity per stirring position (H₂O): 0.4 l

Distance between stirring places: 90 mm

Speed range: 0 to 1000 rpm

Heating rate (1l H₂O in H15): 3 K/min

Temperature range heating plate: RT to 120 °C

Max. temperature medium (dep. on vessel): 70 °C

Temperature display: LED

Heat control accuracy: ±1 K

Speed display: LED line

Set-up plate material: aluminium alloy

Protection class acc. to DIN EN 60529: IP 40

Power supply: 100-240 V, 50/60 Hz



Code	Description	Packaging
LLG09720607	Multi-position magnetic stirrer RT 10 S 2 digital, with 10 stirring places, with heating, with UK plug	1 pz.
LLG09720608	Multi-position magnetic stirrer RT 10 S 3 digital, with 10 stirring places, with heating, with Swiss plug	1 pz.
LLG09720603	Multi-position magnetic stirrer RT 5 digital, with 5 stirring places, with heating	1 pz
LLG09720604	Multi-position magnetic stirrer RT 5 S 2 digital, with 5 stirring places, with heating, with UK plug	1 pz.
LLG09720605	Multi-position magnetic stirrer RT 5 S 3 digital, with 5 stirring places, with heating, with Swiss plug	1 pz.
LLG09720611	Multi-position magnetic stirrer RT 15 S 3 digital, with 15 stirring palces, with heating, with Swiss plug	1 pz.
LLG09720609	Multi-position magnetic stirrer RT 15 digital, with 15 stirring places, with heating	1 pz
LLG09720610	Multi-position magnetic stirrer RT 15 S 2 digital, with 15 stirring places, with heating, with UK plug	1 pz.
LLG09720606	Multi-position magnetic stirrer RT 10 digital, with 10 stirring places, with heating	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

