



LLG-Laboratory coat, 100% cotton

LLG-Laboratory coat, 100% cotton

- Very solid material, excellent wearing comfort
- Snap fasteners (press buttons); easy and quick to open
- With fold-down collar
- One outside breast pocket and two side pockets
- Washable up to 60°C
- Ladies coat with stitched-down back belt

Material: 100% cotton

Colour: White

Weight: 220 g/m²



Code	Description	Packaging
LLG09414346	Laboratory coat, 100% cotton, Mens coat, Size 48	1 pz.
LLG09414343	Laboratory coat, 100% cotton, Ladies coat, Size 40/42	1 pz.
LLG09414344	Laboratory coat, 100% cotton, Ladies coat, Size 44/46	1 pz.
LLG09414345	Laboratory coat, 100% cotton, Mens coat, Size 44/46	1 pz.
LLG09414347	Laboratory coat, 100% cotton, Mens coat, Size 50	1 pz.
LLG09414348	Laboratory coat, 100% cotton, Description Mens coat , Size 52	1 pz.
LLG09414349	Laboratory coat, 100% cotton, Description Mens coat , Size 54	1 pz.
LLG09414350	Laboratory coat, 100% cotton, Mens coat, Size 56	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

