

**Tube roller Roller 6 basic / 6 digital / 10 basic / 10 digital****IKA****Tube roller Roller 6 basic / 6 digital / 10 basic / 10 digital**

Roller shaker with 6 or 10 rolls, with smooth rocking and rolling action and a fixed speed. Particularly for mixing blood samples, solid and liquid suspensions and viscous samples. The unit can be used for all conventional tubes and cylindrical bottles. The easy removable rolls allow simple and quick cleaning of the unit from spilled samples.

- Usable from 4 to 50°C in incubators
- Suitable for continuous operation
- Removable rolls make it suitable for larger tubes
- Easy operation by foil keypad
- Easy removable rolls

Roller 6 digital / 10 digital only:

- Timer: Countdown, adjustable from 1s to 99h 59min
- Counter: Display of shaking time
- Error code display

Specifications Roller 6 /10 basic // Roller 6 /10 digital**Type of movement: rocking and rolling****Amplitude: 24.5 (height)****Speed:****Number of rollers: 30 rpm // 5 ... 80 rpm
6 // 10****Load max.: 2 kg****Timer: - // 1 s ... 99 h 59 min****Ambient temperature: +4°C to +50°C****Max. humidity: 80%****Dimensions (W x D x H):****Roller size (l x Ø): 240 x 545 x 115 mm // 380 x 545 x 115 mm
310 x 32 mm****Weight: 4 kg // 6 kg****Power supply: 100-240V, 50/60Hz****IP code: IP 21**

Code	Description	Packaging
LLG09729970	Roller 6 basic 6 rolls, 30 rpm	1 pz.
LLG09729972	Roller 10 basic 10 rolls, 30 rpm	1 pz.
LLG09729971	Roller 6 digital 6 rolls, 0-80 rpm with timer display	1 pz.
LLG09729973	Roller 10 digital 10 rolls, 0-80 rpm with timer display	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

