



LLG-Desiccators, soda-lime glass, with lid, stopcock and porcelain plate

LLG Labware

LLG-Desiccators, soda-lime glass, with lid, stopcock and porcelain plate

Desiccator made of soda-lime glass with lid, porcelain plate and stopcock. The stopcock has an outer diameter of about 9 mm. Not suitable for vacuum-tight applications.



Code	Description	Packaging
LLG09042740	Desiccators, soda-lime glass, with lid, stopcock and porcelain plate, Ext. diam. 215 mm, Plate diam. 140 mm, Height 224 mm	1 pz
LLG09042741	Desiccators, soda-lime glass, with lid, stopcock and porcelain plate, Ext. diam. 262 mm, Plate diam. 190 mm, Height 239 mm	1 pz
LLG09042742	Desiccators, soda-lime glass, with lid, stopcock and porcelain plate, Ext. diam. 293 mm, Plate diam. 235 mm, Height 278 mm	1 pz
LLG09042743	Desiccators, soda-lime glass, with lid, stopcock and porcelain plate, Ext. diam. 380 mm, Plate diam. 280 mm, Height 414 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

