



Cuvettes for photometers

Lovibond

Cuvettes for photometers



Code	Description	Packaging
LLG09699240	Cuvettes for photometers, Type Set of 10 round cuvettes with cap, Diam. 16 mm, 16mm dia.	10 pz.
LLG06230635	Cuvettes for photometers, Type Set of 5 round cuvettes with cap, Diam. 24 mm, 24 mm dia.	5 pz.
LLG09699231	Cuvettes for photometers, Type Set of 12 round cuvettes with cap, Diam. 24 mm, 24mm dia.	12 pz.
LLG09699243	Cuvettes for photometers, Type Sealing ring for cuvette, Diam. mm, 24mm dia.	12 pz.
LLG09699228	Cuvettes for photometers, Type Adapter, grey, Diam. mm, for 13mm dia. cuvettes	1 pz.
LLG09699226	Cuvettes for photometers, Type Adapter, grey, Diam. mm, for 16mm dia. cuvettes	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

