



Concentrates for ultrasonic baths TICKOPUR

Concentrates for ultrasonic baths TICKOPUR

Concentrates for SONOREX ultrasonic cleaning baths. Enhances cavitation, cleans intensively, kind to most materials and free from chloride. For use cold or warmed up.

TICKOPUR TR 3: Special cleaner, concentrate. Properties: weak acid pH 3.0 at 1%, dosage 5%. Based on: citric acid, surfactants, without phosphate. For: steel, stainless steel, non-ferrous, precious and light metals, glass, ceramics, plastics, rubber. Removes: mineral residues, flash rust, grease, oils, waxes, pigments, grinding, polishing and lapping residues.

TICKOPUR TR 13*: Intensive cleaner, concentrate. Properties: alkaline pH 11.9 at 1 %, dosage 0.1 to 10 %. Based on: surfactants, alkalis. For: steel, stainless steel, precious metals, glass, ceramics, plastics, rubber. Not for: light metals, tin and zinc. Against: Resins, coking residues, soot, greases, oils, waxes, pigments, colour fog, grinding, polishing and lapping residues.

TICKOPUR R 27*: Special cleaner, concentrate. Properties: acid pH 1.9 at 1% solution, dosage 5%. Based on: phosphoric acid. For: steel, stainless steel, precious metal, glass, ceramic rubber, plastic. Removes: rust, lime scale, mineral residues, silicate, phosphate, cement etc.

TICKOPUR R 30: Neutral cleaner, concentrate. Properties: neutral, pH 7 dosage 1 to 5%. Based on: surfactants. For: metal, light metal, glass, ceramics, plastic, rubber. Against: drilling, grinding, polishing and lapping residues, oily and greasy soiling and soot.

TICKOPUR R 33: Universal cleaner, concentrate. Properties: mildly alkaline pH 9.9 at 1%, dosage 3 to 5%. For: metal, light metals, glass, ceramics, plastic, rubber, e-filters, respirators. Against: grinding, polishing and lapping residues, oily and greasy soiling, waxes, silicone films, ink.

TICKOPUR RW 77: Special cleaner with ammonia, concentrate. Properties: mildly alkaline, pH 9.9 at 1%, dosage 5%. Based



on: surfactants, ammonia. For: metal, glass, ceramics, plastic, rubber, boards, watches, test sieves. Beware of light metals. Against: resin, soot, pigments, polishing and lapping residues, greases, oils, waxes, stains, flux.

TICKOPUR J 80 U: Deoxidiser, ready-to-use. Properties: weakly acidic, pH 1.5. Only for use in the application vessel. For: precious metals and their alloys, copper, bronze, brass. Not for light metals, steel, tin, zinc.

Code	Description	Packaging
LLG09192601	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 33, 1 Litres	1 pz.
LLG09192611	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR RW 77, 1 Litres	1 pz.
LLG09192641	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 30, 1 Litres	1 pz.
LLG09192751	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR J 80 U, 1 Litres	1 pz.



Code	Description	Packaging
LLG06301466	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR TR 3, 2 Litres	1 pz.
LLG09192602	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 33, 2 Litres	1 pz.
LLG09192612	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR RW 77, 2 Litres	1 pz.
LLG09192646	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 30, 2 Litres	1 pz.
LLG09192711	Concentrates for ultrasonic baths STAMMOPUR, Type STAMMOPUR 24*, 2 Litres	1 pz.
LLG09192752	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR J 80 U, 2 Litres	1 pz.
LLG09192605	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 33, 5 Litres	1 pz.
LLG09192615	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR RW 77, 5 Litres	1 pz.
LLG09192642	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 30, 5 Litres	1 pz.
LLG09192712	Concentrates for ultrasonic baths STAMMOPUR, Type STAMMOPUR 24*, 5 Litres	1 pz.
LLG09192755	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR J 80 U, 5 Litres	1 pz.
LLG09192761	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 27*, 1 Litres	1 pz.
LLG09192762	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 27*, 2 Litres	1 pz.
LLG09192765	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR R 27*, 5 Litres	1 pz.
LLG09192672	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR TR 13*, 5 Litres	1 pz.
LLG06802927	Concentrates for ultrasonic baths TICKOPUR, Type TICKOPUR TR 3, 1 Litres	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

