

**LLG-Wide-neck jars, glass****LLG Labware****LLG-Wide-neck jars, glass**

Twist-off jar, with wide-neck. Suitable for storing or transporting liquids and coarse solids. The large opening facilitates cleaning and filling of the clear glass.

Please order caps separately.



Code	Description	Packaging
LLG09070301	Wide-neck jars, glass, Nominal capacity 106 ml	12 pz
LLG09070303	Wide-neck jars, glass, Nominal capacity 375 ml	12 pz
LLG09070304	Wide-neck jars, glass, Nominal capacity 720 ml	12 pz
LLG09070305	Wide-neck jars, glass, Nominal capacity 1062 ml	12 pz
LLG09070302	Wide-neck jars, glass, Nominal capacity 212 ml	12 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

