


**Caps, metal, for LLG-Wide-Neck jars****LLG Labware****Caps, metal, for LLG-Wide-Neck jars**

White. Twist-off screw cap for LLG-Wide-neck jars



Code	Description	Packaging
LLG09070306	Caps, metal, for Wide-Neck jars, For jars 106 ml, Sterile -, Lid Metal, Colour white, Diam. 48 mm	12 pz
LLG09070307	Caps, metal, for Wide-Neck jars, For jars 212 / 375 ml, Sterile -, Lid Metal, Colour white, Diam. 63 mm	12 pz
LLG09070308	Caps, metal, for Wide-Neck jars, For jars 720 / 1062 ml, Sterile -, Lid Metal, Colour white, Diam. 82 mm	12 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

