

**pH electrode InLab Semi-Micro**

Mettler-Toledo Online

pH electrode InLab® Semi-Micro

pH-combination electrode with glass shaft, ARGENTALä reference system and S7 electrode head.

Cable - please order separately.

Specifications

Temperature range: 0 ... 100 °C

pH range: 0 ... 12 pH

Temperature sensor: no

Length / Diam.: 100 mm / 6 mm

Electrolyte: XEROLYT®EXTRA Polymer

Diaphragm: open junction

Membrane resistance at 25 °C: 600 MW



Code	Description	Packaging
LLG09041565	pH electrode InLab® Semi-Micro, Type InLab® Semi-Micro	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

