

**pH electrode InLab Micro Pro ISM**

Mettler-Toledo Online

pH electrode InLab® Micro Pro ISM

pH-combination electrode with glass shaft, ARGENTALä reference system with Ag+ trap and MultiPinä electrode head.

Cable - please order separately.

Specifications

Temperature range: 0 ... 100 °C

pH range: 0 ... 14 pH

Temperature sensor: NTC 30 kW

Length / Diam.: 130 mm / 5 mm

Electrolyte: Fluid, KCl 3mol/l

Diaphragm: Ceramic

Membrane resistance at 25 °C: 300 MW



Code	Description	Packaging
LLG06243633	pH electrode InLab® Micro Pro ISM, Type InLab® Micro Pro ISM	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

