

**pH electrode InLab Surface**

Mettler-Toledo Online

pH electrode InLab® Surface

pH-combination electrode with polysulfone shaft, Ag/AgCl reference system and S7 electrode head.

Cable - please order separately.

Specifications

Temperature range: 0 ... 50 °C

pH range: 1 ... 11 pH

Temperature sensor: no

Length / Diam.: 120 mm / 12 mm

Electrolyte: Fluid, KCl 3mol/l, AgCl saturated

Diaphragm: Safelockä ABS ring

Membrane resistance at 25 °C: 800MW



Code	Description	Packaging
LLG09041557	pH electrode InLab® Surface, Type InLab® Surface	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

