

**pH electrode InLab Science Pro ISM**

Mettler-Toledo Online

pH electrode InLab® Science Pro ISM

pH-combination electrode with glass shaft, ARGENTALä reference system with Ag⁺ trap and MultiPinä electrode head.

Cable - please order separately.

Specifications

Temperature range: 0 ... 100 °C

pH range: 0 ... 12 pH

Temperature sensor: NTC 30 kW

Length / Diam.: 170mm / 12mm

Electrolyte: Fluid, KCl 3mol/l

Diaphragm: movable glass sleeve

Membrane resistance at 25 °C: 600 MW



Code	Description	Packaging
LLG09041472	pH electrode InLab® Science Pro ISM, Type InLab® Science Pro ISM	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

