



Precision ground thermometer, enclosed-scale type

Precision ground thermometer, enclosed-scale type

Acc. to DIN 12784, with standard ground joint. Enclosed-scale type, capillary: **transparent prismatic**, diam. top: 11 mm (± 0.5 mm), suitable for calibration.



Code	Description	Packaging
LLG09236601	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236605	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236602	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236606	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236610	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236603	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236607	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236611	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236608	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236613	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236614	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236604	Precision ground thermometer, enclosed-scale type	1 pz.
LLG09236612	Precision ground thermometer, enclosed-scale type	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



