



Sample bags Whirl-Pak, PE, sterile

Sample bags Whirl-Pak®, PE, sterile

Whirl-Pak® bags are sterile polyethylene laboratory sample bags used to transport samples for biological testing and other applications. Available with or without write-on label. With integrated wire closing system. Temperature range: -210 to 82 °C.



Code	Description	Packaging
LLG06091055	Sample bags Whirl-Pak 58 ml, dim. (D x W) 75 x 125 mm, Thickness 0,057 mm, with write-on label	500 pz.
LLG07058005	Sample bags Whirl-Pak 58 ml, dim. (D x W) 75 x 125 mm, Thickness 0,057 mm, with write-on label	500 pz
LLG06266305	Sample bags Whirl-Pak 3637 ml, dim. (D x W) 254 x 508 mm, Thickness 0,102 mm, with write-on label	250 pz.
LLG07981345	Sample bags Whirl-Pak 207 ml, dim. (D x W) 95 x 180 mm, Thickness 0,076 mm, without write-on label	500 pz.
LLG07077212	Sample bags Whirl-Pak 118 ml, dim. (D x W) 75 x 185 mm, Thickness 0,057 mm, with write-on label	500 pz
LLG06305117	Sample bags Whirl-Pak 2041 ml, dim. (D x W) 190 x 380 mm, Thickness 0,076 mm, without write-on label	500 pz
LLG06242083	Sample bags Whirl-Pak 710 ml, dim. (D x W) 150 x 230 mm, Thickness 0,076 mm, with write-on label	500 pz.
LLG06266304	Sample bags Whirl-Pak 2721 ml, dim. (D x W) 254 x 380 mm, Thickness 0,102 mm, with write-on label	250 pz.
LLG07058120	Sample bags Whirl-Pak 532 ml, dim. (D x W) 115 x 230 mm, Thickness 0,064 mm, with write-on label	500 pz
LLG07080208	Sample bags Whirl-Pak 29 ml, dim. (D x W) 65 x 125 mm, Thickness 0,057 mm, with write-on label	500 pz
LLG06266302	Sample bags Whirl-Pak 1242 ml, dim. (D x W) 150 x 380 mm, Thickness 0,076 mm, without write-on label	500 pz.
LLG06280246	Sample bags Whirl-Pak 798 ml, dim. (D x W) 127 x 305 mm, Thickness 0,076 mm, without write-on label	
LLG09404234	Sample bags Whirl-Pak 118 ml, dim. (D x W) 75 x 185 mm, Thickness 0,057 mm, without write-on label	500 pz
LLG09404235	Sample bags Whirl-Pak 532 ml, dim. (D x W) 115 x 230 mm, Thickness 0,064 mm, without write-on label	500 pz
LLG09404236	Sample bags Whirl-Pak 1065 ml, dim. (D x W) 125 x 380 mm, Thickness 0,076 mm, without write-on label	500 pz





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

