

**Nessler Tube, without spout, borosilicate 3.3****NORMAX****Nessler Tube, without spout, borosilicate 3.3**

For use, for example, with the Nessleriser systems from Lovibond®.



Code	Description	Packaging
LLG09170622	Nessler Tube, without spout, borosilicate 3.3, tall form, graduated at 100 ml, 100 ml, dim. (dia. x H) 24 x 375 mm	1 pz
LLG09170623	Nessler Tube, without spout, borosilicate 3.3, short form, graduated at 50 ml, 50 ml, dim. (dia. x H) 25 x 175 mm	1 pz
LLG09170624	Nessler Tube, without spout, borosilicate 3.3, short form, graduated at 100 ml, 100 ml, dim. (dia. x H) 32 x 200 mm	1 pz
LLG09170621	Nessler Tube, without spout, borosilicate 3.3, tall form, graduated at 50 ml, 50 ml, dim. (dia. x H) 20 x 300 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

