

**Sample bags Whirl-Pak Stand-Up, PE, sterile, free standing****Sample bags Whirl-Pak® Stand-Up, PE, sterile, free standing**

Easy sample handling due to the specially gusseted bottom which flattens out when the sample is placed inside, allowing the bag to stand up on its own. It works just like a bottle or beaker, however it cost less, is shatterproof and uses a fraction of the storage space. The integrated wire closing system ensures the tightness of the sample bag. With write-on label for clean and easy identification. Temperature range: -210 to 82 °C.



Code	Description	Packaging
LLG06260542	Sample bags Whirl-Pak® Stand-Up, PE, sterile, free standing, Nominal capacity 1065 ml, dim. (D x W) 125 x 380 mm, Thickness 0,102 mm	250 pz.
LLG06266303	Sample bags Whirl-Pak® Stand-Up, PE, sterile, free standing, Nominal capacity 710 ml, dim. (D x W) 150 x 230 mm, Thickness 0,076 mm	500 pz.
LLG06260543	Sample bags Whirl-Pak® Stand-Up, PE, sterile, free standing, Nominal capacity 2041 ml, dim. (D x W) 190 x 380 mm, Thickness 0,102 mm	250 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

