


**Laboratory jack LABLIFT m****IKA****Laboratory jack LABLIFT m**

IKA® lab lifts are suitable for raising and supporting small laboratory devices or vessels. The auxiliary springs ensures that users can raise large loads easily, especially from the lowest position.

- The hand wheel ensures that the operation is smooth and comfortable
- The self-lubricating bearing significantly reduces friction



Code	Description	Packaging
LLG09720419	Laboratory jack LABLIFT m, Width 200 mm, Depth 310 mm, Min. height 70 mm, Max. height 300 mm	1 pz



CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System

