

**Laboratory jack LABLIFT m****IKA****Laboratory jack LABLIFT m**

IKA® lab lifts are suitable for raising and supporting small laboratory devices or vessels. The auxiliary springs ensures that users can raise large loads easily, especially from the lowest position.

- The hand wheel ensures that the operation is smooth and comfortable
- The self-lubricating bearing significantly reduces friction



Code	Description	Packaging
LLG09720419	Laboratory jack LABLIFT m, Width 200 mm, Depth 310 mm, Min. height 70 mm, Max. height 300 mm	1 pz



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
Reagents operates with
a Certified Quality
Management System

