



**Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C**

**Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C**

**Exchangeable blocks:**

Flexibility has never been this easy. Eppendorf offers a variety of SmartBlocks for tubes from 0.5ml to 50ml.

- Fast and simple block exchange due to Eppendorf Quick Release
- Optimized block design guarantees maximum temperature transfer to the sample

**Transfer Racks:**

- Autoclavable transport aid for tubes from 0.5 to 2.0 ml
- For transferring up to 24 tubes in / from the heating block at the same time
- Offers protection against burns when transferring hot samples

**SmartExtender:**

A heatable frame/attachment for mixer, with separate temperature control (3 °C above RT to 110 °C), for simultaneous operation at two different temperatures.

- Enables 2 incubations at the same time



| Code        | Description  | Packaging |
|-------------|--|-----------|
| LLG06280936 | Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C, Type Transfer Rack, For 24 tubes, 0.5 ml         |           |
| LLG06280935 | Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C, Type Transfer Rack, For 24 tubes, 1.5 ... 2.0 ml |           |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

