



Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C

Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C

Exchangeable blocks:

Flexibility has never been this easy. Eppendorf offers a variety of SmartBlocks for tubes from 0.5ml to 50ml.

- Fast and simple block exchange due to Eppendorf Quick Release
- Optimized block design guarantees maximum temperature transfer to the sample

Transfer Racks:

- Autoclavable transport aid for tubes from 0.5 to 2.0 ml
- For transferring up to 24 tubes in / from the heating block at the same time
- Offers protection against burns when transferring hot samples

SmartExtender:

A heatable frame/attachment for mixer, with separate temperature control (3 °C above RT to 110 °C), for simultaneous operation at two different temperatures.

- Enables 2 incubations at the same time



Code	Description	Packaging
LLG06280936	Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C, Type Transfer Rack, For 24 tubes, 0.5 ml	
LLG06280935	Exchangeable blocks Eppendorf SmartBlocks and accessories for Eppendorf ThermoMixer C and ThermoStat C, Type Transfer Rack, For 24 tubes, 1.5 ... 2.0 ml	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

