



Test tube racks, PP

Test tube racks, PP

Blue. For tubes up to 18mm in diameter. Autoclavable and unbreakable. Can be dismantled. Can be recycled.



Code	Description	Packaging
LLG09193112	Test tube racks, PP, For tubes diam. 18 mm, Array 12	1 pz
LLG09193124	Test tube racks, upto 18 mm diam, 24 tube	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

