




Test tube racks, wire/nylon

Test tube racks, wire/nylon

Wire with white nylon coating. 3-tier design, second tier approximately 30 mm high. Opening size 18 x 18 mm. Temperature resistant up to 80 °C max.



Code	Description	Packaging
LLG09193776	Test tube racks, wire/nylon, Height 70 mm, Array 10 x 10	1 pz.
LLG09193764	Test tube racks, wire/nylon, Height 70 mm, Array 2 x 10	1 pz.
LLG09193766	Test tube racks, wire/nylon, Height 70 mm, Array 2 x 12	1 pz.
LLG09193762	Test tube racks, wire/nylon, Height 70 mm, Array 2 x 6	1 pz.
LLG09193768	Test tube racks, wire/nylon, Height 70 mm, Array 3 x 12	1 pz.
LLG09193774	Test tube racks, wire/nylon, Height 70 mm, Array 4 x 12	1 pz.
LLG09193788	Test tube racks, wire/nylon, Height 100 mm, Array 4 x 12	1 pz.
LLG09193772	Test tube racks, wire/nylon, Height 70 mm, Array 4 x 6	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

