


**Wire gauzes****RSG****Wire gauzes**

Plain, stainless steel wire with heat resistant center, mesh: 0.6mm x 2.0mm



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09033653 | Wire gauzes, Width 135 mm, Depth 135 mm | 1 pz.     |
| LLG09033655 | Wire gauzes, Width 155 mm, Depth 155 mm | 1 pz.     |
| LLG09033657 | Wire gauzes, Width 175 mm, Depth 175 mm | 1 pz.     |
| LLG09033659 | Wire gauzes, Width 195 mm, Depth 195 mm | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

