

**Portable Desiccator DURAPorter, PC****Heathrow Scientific****Portable Desiccator DURAPorterä, PC**

Ideal for transporting biological and clinical samples, as well as instruments and products that you want to keep clean and dry under testing environmental conditions. Are light weight and easy to carry; they are ideal for moving samples that must be protected from humidity, dust and other environmental influence.

The clarity of the DURAPorterä enables the user to visibly check sample integrity, so that it can be opened using the right precautions.

- Polycarbonate construction are available in three colours for use/owner identification
- Silicone water-tight seal and three lid clasps provide secure closure
- Carry handle folds neatly into body space allowing the units to be stacked
- Lid opens a full 180°, allowing total access to contents and making it easier to clean
- Will accommodate 13 and 16 mm tubes in a 72-place one rack®, or similar sized/brand of rack
- Use the separators provided to create your own compartment sizes, or use empty to transport larger products and equipment
- Autoclavable

**Matching rack (e.g. Cat. No. 9193 997, 9193 998) please order separately.
Further colours available on request.**



Code	Description	Packaging
LLG09194005	Portable Desiccator DURAPorter, PC, Type DURAPorterTM, Colour Clear with blue handles, dim. (l x w x h) 380 x 196 x 160 mm, Material PC	1 pz.
LLG09194007	Portable Desiccator DURAPorter, PC, Type DURAPorterTM, Colour Red with red handles, dim. (l x w x h) 380 x 196 x 160 mm, Material PC	1 pz.
LLG09194006	Portable Desiccator DURAPorter, PC, Type DURAPorterTM, Colour Yellow with yellow handles, dim. (l x w x h) 380 x 196 x 160 mm, Material PC	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

