

**Inoculating needles smartPicks, wood****Heathrow Scientific****Inoculating needles smartPicks®, wood**

Available in three lengths and two needle sizes to suit a wide range of application needs. The Mini fits neatly into a 1.5 mL microfuge tube allowing the lid to be closed for full sample transfer/recovery.

- Naturally biodegradable White Birch wood
- Manufactured in an FDA-registered factory and complies with FDA standards for food contact
- Autoclavable and available in semi or needle point styles for delicate operations



Code	Description	Packaging
LLG06257177	Inoculating needles smartPicks®, wood, Type smartPicks® Mini, Needle point, for 1.5ml tubes	2000 pz.
LLG06257176	Inoculating needles smartPicks®, wood, Type smartPicks® Mini, Semi point, for 1.5ml tubes	2000 pz.
LLG06257178	Inoculating needles smartPicks®, wood, Type smartPicks® Standard, Semi point, single tip	1500 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

