

**Tube Rack Origami, PP****Tenak****Tube Rack Origami, PP**

The Origami Tube Rack is designed to be height adjustable which makes it very flexible and usable for almost any test tube on the market. Depending on the hole diameter the racks can be used for microcentrifuge tubes, cryovials, vacutainers, and 15ml and 50ml falcon tubes. The Origami Tube Rack has a printed alphanumeric grid for easy location of the samples. It is made of polypropylene and can be used in temperatures from -45 °C to 60 °C. Best choice for mobile Labs. After usage demount it and store it in a drawer.



Code	Description	Packaging
LLG09194017	Tube Rack Origami, PP, For tubes diam. 30 mm, Array 15, Dimensions (W x D x H) 240 x 120 x 140 mm	5 pz.
LLG09194015	Tube Rack Origami, PP, For tubes diam. 12.5 mm, Array 45, Dimensions (W x D x H) 240 x 120 x 140 mm	5 pz.
LLG09194016	Tube Rack Origami, PP, For tubes diam. 18 mm, Array 28, Dimensions (W x D x H) 240 x 120 x 140 mm	5 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

