

**Distilling Adapters, borosilicate glass 3.3****ISOLAB****Distilling Adapters, borosilicate glass 3.3**

Manufactured from borosilicate glass 3.3. Suitable to be used in distillation assemblies for connecting flasks to condenser. Upper socket designed to accept thermometer either by screw thread adapter or by thermometer socket.

With side cone at 75° angled to the center socket.



Code	Description	Packaging
LLG09043301	Distilling adaptersocket NS 14/23, cone	1 pz
LLG09043305	Distilling adaptersocket NS 14/23, cone	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

