

**LLG-Labware pipette tips, economy 2.0, PP****LLG Labware****LLG-Labware pipette tips, economy 2.0, PP**

The new generation of LLG Labware Pipette tips economy.

Available in bulk quantities in bags or in easy-to-open boxes with color coded racks. Racked-tip boxes are delivered in a dispenser, which can be opened from two sides facilitating the use of the product on two work spaces. Compatible with all common pipette brands. You can find the appropriate compatibility list on [www.llg-labware.com](http://www.llg-labware.com) in the Downloads section.

- Improved optical clarity and pipette fit
- Colored coded rack inserts for easy identification
- Graduated (5 and 10 ml tips without graduation)
- Autoclavable
- DNase and RNase free



Code	Description	Packaging
LLG04668775	LLG TIPS ECO2 0.1-10 $\mu$ L clear bulk 1000	1000 pz
LLG04668776	LLG TIPS ECO2 0.1-10 $\mu$ L clear rack 10x96	960 pz
LLG04668778	LLG TIPS ECO2 1-200 $\mu$ L clear rack 10x96	10x96
LLG04668779	LLG TIPS ECO2 1-200 $\mu$ L yellow rack 10x96	960 pz
LLG04668780	LLG TIPS ECO2 100-1000 $\mu$ L clear bulk 1000	1000 pz
LLG04668782	LLG TIPS ECO2 100-1000 $\mu$ L clear rack 8x96	768 pz
LLG04668783	LLG TIPS ECO2 100-1000 $\mu$ L blue rack 8x96	768 pz
LLG09409097	LLG-Pipette tips Economy 100-5000 $\mu$ /250	250 pz
LLG09409098	LLG-Pipette tips Economy 1000-10000/100	100 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

