



Receiving flasks, DURAN

Receiving flasks, DURAN®

Especially for rotary evaporators, with spherical joint socket S 35, with and without transparent PUR coating, autoclavable, material: DURAN®.



Code	Description	Packaging
LLG09196224	Receiving flasks, DURAN®, 100 ml, socket 35 NS, Flask diam. 64 mm, uncoated	1 pz.
LLG09196254	Receiving flasks, DURAN®, 1000 ml, socket 35 NS, Flask diam. 131 mm, uncoated	1 pz.
LLG06200782	Receiving flasks, DURAN®, 2000 ml, socket 35 NS, Flask diam. 166 mm, PUR coating	1 pz.
LLG07653714	Receiving flasks, DURAN®, 500 ml, socket 35 NS, Flask diam. 105 mm, PUR coating	1 pz.
LLG09812217	Receiving flasks, DURAN®, 1000 ml, socket 35 NS, Flask diam. 131 mm, PUR coating	1 pz.
LLG09196236	Receiving flasks, DURAN®, 250 ml, socket 35 NS, Flask diam. 85 mm, uncoated	1 pz.
LLG09196244	Receiving flasks, DURAN®, 500 ml, socket 35 NS, Flask diam. 105 mm, uncoated	1 pz.
LLG09196263	Receiving flasks, DURAN®, 2000 ml, socket 35 NS, Flask diam. 166 mm, uncoated	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

