

**Disposable tweezers, PS, blue****Bürkle GmbH****Disposable tweezers, PS, blue**

For precise laboratory work, for gripping and holding small objects. The blue products, or parts of them, can be quickly found and easily identified in a visual check, if they ever get into production.

- Serrated handles
- Sterile
- Cleanroom manufactured (cleanroom class 7)
- Food-safe

Form Pointed / angled: fine, smooth, curved gripping tip

Length 130 mm

Sterile +



Code	Description	Packaging
LLG04670767	Disposable tweezers, PS, blue, Form Wide tip: smooth and blunt tip , Length 130 mm, Sterile +	100 pc
LLG04670765	Disposable tweezers, PS, blue, Form Pointed / angled: fine, smooth, curved gripping tip , Length 130 mm, Sterile +	100 pc
LLG04670766	Disposable tweezers, PS, blue, Form Pointed: straight, inside, serrated tip , Length 130 mm, Sterile +	100 pc



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

