

**OTAL disposable pumps with PE conveying pipes****Bürkle GmbH****OTAL® disposable pumps with PE conveying pipes**

Space-saving container pump for single use.

The medium to be filled only comes into contact with the conveying pipe.

The conveying pipe is made of one piece and can easily be shortened to the respective container height with the supplied cutter for good residual emptying.

The disposable conveying pipe can be replaced quickly and easily after use.

Scope of supply: Pump ball, universal plug made of soft PVC for container openings of 30 - 70 mm Ø, cutter and 10 disposable conveying pipes



Code	Description	Packaging
LLG04670768	OTAL® disposable foot pump , Ø 12 mm, Immersion depth 90 mm, Flow rate 14 l / min.	
LLG04670769	OTAL®disposable hand pump , Ø 12 mm, Immersion depth 90 mm, Flow rate 14 l / min.	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

